

Product Datasheet

Human LIPD ELISA Kit (orb440552)

Description	Human Lipase, Lipoprotein (LIPD) Sandwich ELISA Kit
Reactivity	Human
Range	23.437-1500pg/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma, tissue homogenates, cell lysates and other biological fluids.
Sample Types	serum, plasma, tissue homogenates, cell lysates and other biological fluids.
Uniprot ID	P06858
Sensitivity	8.65pg/mL
Entrez	4023
Expiration Date	Please enquire.