


Product Datasheet

HLA-A antibody (orb44056)

Description	Mouse monoclonal antibody to HLA-A
Reactivity	Bovine, Human, Primate
Conjugation	Unconjugated
Tested Applications	ELISA, FC, IP
Immunogen	COS-7 African green monkey kidney cells
Target	HLA Class I
Preservatives	Tris buffered saline (TBS), pH 8.0, 15 mM sodium azide
Concentration	1 mg/ml
Storage	Store at 2-8°C. Do not freeze.
Note	For research use only
Application notes	Flow cytometry: Recommended dilution: 4 µg/ml.
Isotype	Mouse IgG3
Clonality	Monoclonal
Purity	Purified by protein-A affinity chromatography.
Dilution Range	Flow cytometry: Recommended dilution: 4 µg/ml.
Expiration Date	12 months from date of receipt.



Flow
cytometry
analysis
(surface
stainin...