

Product Datasheet

Mouse LOXL1 ELISA Kit (orb440583)

Description	Mouse Lysyl Oxidase Like Protein 1 (LOXL1) Sandwich ELISA Kit
Reactivity	Mouse
Range	31.25-2000pg/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma, tissue homogenates and other biological fluids.
Sample Types	serum, plasma, tissue homogenates and other biological fluids.
Uniprot ID	P97873
Sensitivity	14.8pg/mL
Entrez	16949
Expiration Date	Please enquire.