

Product Datasheet

Human LPCAT2 ELISA Kit (orb440597)

Description	Human Lysophosphatidylcholine Acyltransferase 2 (LPCAT2) Sandwich ELISA Kit
Reactivity	Human
Range	1.56-100ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma, tissue homogenates, cell lysates, cell culture supernates and other biological fluids.
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernatants and other biological fluids.
Uniprot ID	Q7L5N7
Sensitivity	0.58ng/mL
Entrez	54947
Expiration Date	Please enquire.