

Product Datasheet

Rat LpPLA2 ELISA Kit (orb440610)

Description	Rat Phospholipase A2, Lipoprotein Associated (LpPLA2) Sandwich ELISA Kit
Reactivity	Rat
Range	62.5-4000pg/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma, tissue homogenates and other biological fluids.
Sample Types	serum, plasma, tissue homogenates and other biological fluids.
Uniprot ID	Q5M7T7
Sensitivity	15.8pg/mL
Entrez	301265
Expiration Date	Please enquire.