

Product Datasheet

Rat LRP8 ELISA Kit (orb440623)

Description	Rat Low Density Lipoprotein Receptor Related Protein 8 (LRP8) Sandwich ELISA Kit
Reactivity	Rat
Range	31.25-2000pg/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma and other biological fluids.
Sample Types	serum, plasma and other biological fluids.
Uniprot ID	D3ZE75
Sensitivity	13.3pg/mL
Entrez	362558
Expiration Date	Please enquire.