

Product Datasheet

Human LRRC3C ELISA Kit (orb440624)

Description	Human Leucine Rich Repeat Containing Protein 3C (LRRC3C) Sandwich ELISA Kit
Reactivity	Human
Range	0.156-10ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	tissue homogenates and other biological fluids.
Sample Types	tissue homogenates and other biological fluids.
Uniprot ID	A6NJV4
Sensitivity	0.066ng/mL
Entrez	100505591
Expiration Date	Please enquire.