

Product Datasheet

Mouse LTA4H ELISA Kit (orb440632)

Description	Mouse Leukotriene A4 Hydrolase (LTA4H) Sandwich ELISA Kit
Reactivity	Mouse
Range	0.156-10ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma, tissue homogenates and other biological fluids.
Sample Types	serum, plasma, tissue homogenates and other biological fluids.
Uniprot ID	P24527
Sensitivity	0.051ng/mL
Entrez	16993
Expiration Date	Please enquire.