

Product Datasheet

Human LZM ELISA Kit (orb440645)

Description	Human Lysozyme (LZM) Sandwich ELISA Kit
Reactivity	Human
Range	12.5-800pg/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma and other biological fluids.
Sample Types	serum, plasma and other biological fluids.
Uniprot ID	P61626
Sensitivity	4.8pg/mL
Entrez	4069
Expiration Date	Please enquire.