

Product Datasheet

Human MAN1A1 ELISA Kit (orb440653)

Description	Human Mannosidase Alpha Class 1A Member 1 (MAN1A1) Sandwich ELISA Kit
Reactivity	Human
Range	0.156-10ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	tissue homogenates and other biological fluids.
Sample Types	tissue homogenates and other biological fluids.
Uniprot ID	P33908
Sensitivity	0.054ng/mL
Entrez	4121
Expiration Date	Please enquire.