

Product Datasheet

Human MCL1 ELISA Kit (orb440700)

Description	Human Myeloid Cell Leukemia Sequence 1, Bcl2 Related (MCL1) Sandwich ELISA Kit
Reactivity	Human
Range	0.312-20ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	tissue homogenates, cell lysates and other biological fluids.
Sample Types	tissue homogenates, cell lysates and other biological fluids.
Uniprot ID	Q07820
Sensitivity	0.118ng/mL
Entrez	4170
Expiration Date	Please enquire.