

Product Datasheet

Rabbit MCP1 ELISA Kit (orb440709)

Description	Rabbit Monocyte Chemotactic Protein 1 (MCP1) Sandwich ELISA Kit
Reactivity	Rabbit
Range	31.25-2000pg/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma, tissue homogenates, cell lysates, cell culture supernates and other biological fluids.
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernatants and other biological fluids.
Uniprot ID	P28292
Sensitivity	14.2pg/mL
Entrez	100009130
Expiration Date	Please enquire.