

Product Datasheet

Canine MCP2 ELISA Kit (orb440712)

Description	Canine Monocyte Chemotactic Protein 2 (MCP2) Sandwich ELISA Kit
Reactivity	Canine
Range	31.25-2000pg/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma, tissue homogenates and other biological fluids.
Sample Types	serum, plasma, tissue homogenates and other biological fluids.
Uniprot ID	Q68AY9
Sensitivity	9.7pg/mL
Entrez	448792
Expiration Date	Please enquire.