

Product Datasheet

Mouse MCP2 ELISA Kit (orb440714)

Description	Mouse Monocyte Chemotactic Protein 2 (MCP2) Sandwich ELISA Kit
Reactivity	Mouse
Range	15.625-1000pg/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma, tissue homogenates, cell lysates, cell culture supernates and other biological fluids.
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernatants and other biological fluids.
Uniprot ID	Q9Z121
Sensitivity	6.3pg/mL
Entrez	20307
Expiration Date	Please enquire.