

Product Datasheet

Mouse MECP2 ELISA Kit (orb440736)

Description	Mouse Methyl CpG Binding Protein 2 (MECP2) Sandwich ELISA Kit
Reactivity	Mouse
Range	0.312-20ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	tissue homogenates, cell lysates and other biological fluids.
Sample Types	tissue homogenates, cell lysates and other biological fluids.
Uniprot ID	Q9Z2D6
Sensitivity	0.102ng/mL
Entrez	17257
Expiration Date	Please enquire.