

Product Datasheet

Human MELK ELISA Kit (orb440737)

Description	Human Maternal Embryonic Leucine Zipper Kinase (MELK) Sandwich ELISA Kit
Reactivity	Human
Range	0.312-20ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	tissue homogenates, cell lysates and other biological fluids.
Sample Types	tissue homogenates, cell lysates and other biological fluids.
Uniprot ID	Q14680
Sensitivity	0.105ng/mL
Entrez	9833
Expiration Date	Please enquire.