

Product Datasheet

Mouse MGP ELISA Kit (orb440758)

Description	Mouse Matrix Gla Protein (MGP) Sandwich ELISA Kit
Reactivity	Mouse
Range	78.125-5000pg/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma, tissue homogenates and other biological fluids.
Sample Types	serum, plasma, tissue homogenates and other biological fluids.
Uniprot ID	P19788
Sensitivity	33pg/mL
Entrez	17313
Expiration Date	Please enquire.