

Product Datasheet

Rat MGP ELISA Kit (orb440759)

Description	Rat Matrix Gla Protein (MGP) Sandwich ELISA Kit
Reactivity	Rat
Range	93.75-6000pg/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma, tissue homogenates, cell lysates, cell culture supernates and other biological fluids.
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernatants and other biological fluids.
Uniprot ID	P08494
Sensitivity	39.65pg/mL
Entrez	25333
Expiration Date	Please enquire.