

Product Datasheet

Human MHCDRb1 ELISA Kit (orb440767)

Description	Human Major Histocompatibility Complex Class II DR Beta 1 (MHCDRb1) Sandwich ELISA Kit
Reactivity	Human
Range	1.56-100ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma, tissue homogenates and other biological fluids.
Sample Types	serum, plasma, tissue homogenates and other biological fluids.
Uniprot ID	P01912
Sensitivity	0.72ng/mL
Entrez	3123
Expiration Date	Please enquire.