

Product Datasheet

Mouse MID1 ELISA Kit (orb440773)

Description	Mouse Midline 1 (MID1) Sandwich ELISA Kit
Reactivity	Mouse
Range	0.156-10ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	tissue homogenates and other biological fluids.
Sample Types	tissue homogenates and other biological fluids.
Uniprot ID	070583
Sensitivity	0.052ng/mL
Entrez	17318
Expiration Date	Please enquire.