

Product Datasheet

Human MIOX ELISA Kit (orb440777)

Description	Human Myo Inositol Oxygenase (MIOX) Sandwich ELISA Kit
Reactivity	Human
Range	0.156-10ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	tissue homogenates, cell lysates and other biological fluids.
Sample Types	tissue homogenates, cell lysates and other biological fluids.
Uniprot ID	Q9UGB7
Sensitivity	0.058ng/mL
Entrez	55586
Expiration Date	Please enquire.