

Product Datasheet

Bovine MIP1a ELISA Kit (orb440778)

Description	Bovine Macrophage Inflammatory Protein 1 Alpha (MIP1a) Sandwich ELISA Kit
Reactivity	Bovine
Range	15.625-1000pg/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma, tissue homogenates, cell lysates, cell culture supernates and other biological fluids.
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernatants and other biological fluids.
Uniprot ID	Q8SQA6
Sensitivity	7.4pg/mL
Entrez	282170
Expiration Date	Please enquire.