

Product Datasheet

Canine MIP1a ELISA Kit (orb440779)

Description	Canine Macrophage Inflammatory Protein 1 Alpha (MIP1a) Sandwich ELISA Kit
Reactivity	Canine
Range	78.125-5000pg/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma and other biological fluids.
Sample Types	serum, plasma and other biological fluids.
Uniprot ID	Q68A92
Sensitivity	28pg/mL
Entrez	448787
Expiration Date	Please enquire.