

Product Datasheet

Human Midkine ELISA Kit (orb440804)

Description	Human Midkine (MK) Sandwich ELISA Kit
Reactivity	Human
Range	0.312-20ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	tissue homogenates, cell lysates, cell culture supernates and other biological fluids.
Sample Types	tissue homogenates, cell lysates, cell culture supernatants and other biological fluids.
Uniprot ID	P21741
Sensitivity	0.127ng/mL
Entrez	4192
Expiration Date	Please enquire.