

Product Datasheet

Human MMP10 ELISA Kit (orb440810)

Description	Human Matrix Metalloproteinase 10 (MMP10) Sandwich ELISA Kit
Reactivity	Human
Range	78.125-5000pg/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma and other biological fluids.
Sample Types	serum, plasma and other biological fluids.
Uniprot ID	P09238
Sensitivity	25pg/mL
Entrez	4319
Expiration Date	Please enquire.