

Product Datasheet

Mouse MMP10 ELISA Kit (orb440811)

Description	Mouse Matrix Metalloproteinase 10 (MMP10) Sandwich ELISA Kit
Reactivity	Mouse
Range	6.25-400pg/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma and other biological fluids.
Sample Types	serum, plasma and other biological fluids.
Uniprot ID	O55123
Sensitivity	2.35pg/mL
Entrez	17384
Expiration Date	Please enquire.