

Product Datasheet

Human MMP14 ELISA Kit (orb440829)

Description	Human Matrix Metalloproteinase 14 (MMP14) Sandwich ELISA Kit
Reactivity	Human
Range	1.56-100ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma, tissue homogenates and other biological fluids.
Sample Types	serum, plasma, tissue homogenates and other biological fluids.
Uniprot ID	P50281
Sensitivity	0.63ng/mL
Entrez	4323
Expiration Date	Please enquire.