

Product Datasheet

Mouse MMP19 ELISA Kit (orb440836)

Description	Mouse Matrix Metalloproteinase 19 (MMP19) Sandwich ELISA Kit
Reactivity	Mouse
Range	0.156-10ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	tissue homogenates, cell lysates and other biological fluids.
Sample Types	tissue homogenates, cell lysates and other biological fluids.
Uniprot ID	Q9JH10
Sensitivity	0.054ng/mL
Entrez	58223
Expiration Date	Please enquire.