

Product Datasheet

Porcine MMP2 ELISA Kit (orb440855)

Description	Porcine Matrix Metalloproteinase 2 (MMP2) Sandwich ELISA Kit
Reactivity	Porcine
Range	78.125-5000pg/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma, tissue homogenates, cell lysates and other biological fluids.
Sample Types	serum, plasma, tissue homogenates, cell lysates and other biological fluids.
Uniprot ID	Q95JA4
Sensitivity	33pg/mL
Entrez	397391
Expiration Date	Please enquire.