

Product Datasheet

Porcine MMP3 ELISA Kit (orb440863)

Description	Porcine Matrix Metalloproteinase 3 (MMP3) Sandwich ELISA Kit
Reactivity	Porcine
Range	4.687-300ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma and other biological fluids.
Sample Types	serum, plasma and other biological fluids.
Uniprot ID	C9VZX5
Sensitivity	1.8ng/mL
Entrez	396769
Expiration Date	Please enquire.