

Product Datasheet

Rabbit MMP3 ELISA Kit (orb440865)

Description	Rabbit Matrix Metalloproteinase 3 (MMP3) Sandwich ELISA Kit
Reactivity	Rabbit
Range	0.156-10ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma and other biological fluids.
Sample Types	serum, plasma and other biological fluids.
Uniprot ID	P28863
Sensitivity	0.061ng/mL
Entrez	100009111
Expiration Date	Please enquire.