

Product Datasheet

Bovine MMP7 ELISA Kit (orb440866)

Description	Bovine Matrix Metalloproteinase 7 (MMP7) Sandwich ELISA Kit
Reactivity	Bovine
Range	0.312-20ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma, tissue homogenates and other biological fluids.
Sample Types	serum, plasma, tissue homogenates and other biological fluids.
Uniprot ID	Q148N3
Sensitivity	0.089ng/mL
Entrez	286794
Expiration Date	Please enquire.