

Product Datasheet

Equine MMP9 ELISA Kit (orb440878)

Description	Equine Matrix Metalloproteinase 9 (MMP9) Sandwich ELISA Kit
Reactivity	Equine
Range	15.625-1000pg/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma, tissue homogenates and other biological fluids.
Sample Types	serum, plasma, tissue homogenates and other biological fluids.
Uniprot ID	A8HDR2
Sensitivity	6.7pg/mL
Entrez	100056599
Expiration Date	Please enquire.