

Product Datasheet

Rat MOG ELISA Kit (orb440890)

Description	Rat Myelin Oligodendrocyte Glycoprotein (MOG) Sandwich ELISA Kit
Reactivity	Rat
Range	0.156-10ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	tissue homogenates and other biological fluids.
Sample Types	tissue homogenates and other biological fluids.
Uniprot ID	Q63345
Sensitivity	0.057ng/mL
Entrez	24558
Expiration Date	Please enquire.