

Product Datasheet

Human MPIF2 ELISA Kit (orb440894)

Description	Human Myeloid Progenitor Inhibitory Factor 2 (MPIF2) Sandwich ELISA Kit
Reactivity	Human
Range	32.5-2000pg/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma, tissue homogenates, cell lysates, cell culture supernates and other biological fluids.
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernatants and other biological fluids.
Uniprot ID	000175
Sensitivity	15.2pg/mL
Entrez	6369
Expiration Date	Please enquire.