

Product Datasheet

Equine MPO ELISA Kit (orb440898)

Description	Equine Myeloperoxidase (MPO) Sandwich ELISA Kit
Reactivity	Equine
Range	3.125-200ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma and other biological fluids.
Sample Types	serum, plasma and other biological fluids.
Uniprot ID	A5A4Z8
Sensitivity	1.36ng/mL
Entrez	100070920
Expiration Date	Please enquire.