

Product Datasheet

Human MRC1 ELISA Kit (orb440909)

Description	Human Macrophage mannose receptor 1 (MRC1) Sandwich ELISA Kit
Reactivity	Human
Range	0.312-20ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma, tissue homogenates and other biological fluids.
Sample Types	serum, plasma, tissue homogenates and other biological fluids.
Uniprot ID	P22897
Sensitivity	0.109ng/mL
Entrez	4360
Expiration Date	Please enquire.