

Product Datasheet

Human MSMb ELISA Kit (orb440921)

Description	Human Microseminoprotein Beta (MSMb) Sandwich ELISA Kit
Reactivity	Human
Range	0.156-10ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma, seminal plasma, tissue homogenates, cell lysates, cell culture supernates and other biological fluids.
Sample Types	serum, plasma, seminal plasma, tissue homogenates, cell lysates, cell culture supernatants and other biological fluids.
Uniprot ID	P08118
Sensitivity	0.061ng/mL
Entrez	4477
Expiration Date	Please enquire.