

Product Datasheet

Mouse MT1 ELISA Kit (orb440938)

Description	Mouse Metallothionein 1 (MT1) Sandwich ELISA Kit
Reactivity	Mouse
Range	62.5-4000pg/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma, tissue homogenates and other biological fluids.
Sample Types	serum, plasma, tissue homogenates and other biological fluids.
Uniprot ID	P02802
Sensitivity	35.8pg/mL
Entrez	17748
Expiration Date	Please enquire.