

Product Datasheet

LTBP4 antibody (Biotin) (orb44094)

Description	LTBP4 antibody (Biotin)
Species/Host	Rabbit
Reactivity	Human
Conjugation	Biotin
Tested Applications	ELISA
Immunogen	Recombinant Human Latent-transforming growth factor beta-binding protein 4 protein (977-1180AA)
Target	LTBP4
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	> 95%, Protein G purified
Uniprot ID	Q8N2S1
Dilution Range	ELISA
Expiration Date	12 months from date of receipt.