

Product Datasheet

Mouse MTFMT ELISA Kit (orb440946)

Description	Mouse Mitochondrial Methionyl tRNA Formyltransferase (MTFMT) Sandwich ELISA Kit
Reactivity	Mouse
Range	0.156-10ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	tissue homogenates, cell lysates and other biological fluids.
Sample Types	tissue homogenates, cell lysates and other biological fluids.
Uniprot ID	Q9D799
Sensitivity	0.055ng/mL
Entrez	69606
Expiration Date	Please enquire.