

Product Datasheet

Rat MTL ELISA Kit (orb440949)

Description	Rat Motilin (MTL) Sandwich ELISA Kit
Reactivity	Rat
Range	93.75-6000pg/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma and other biological fluids.
Sample Types	serum, plasma and other biological fluids.
Uniprot ID	A8IRIO
Sensitivity	42.5pg/mL
Entrez	126231
Expiration Date	Please enquire.