


Product Datasheet

CD108 antibody (orb44095)

Description	Mouse Monoclonal to CD108.
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC, IP, WB
Immunogen	HPB-ALL human T cell line
Target	CD108
Preservatives	Tris buffered saline (TBS), pH 8.0, 15 mM sodium azide
Concentration	1 mg/ml
Storage	Store at 2-8°C. Do not freeze.
Note	For research use only
Application notes	Flow cytometry: Recommended dilution: 4 µg/ml. Western blotting: Non-reducing conditions.
Isotype	Mouse IgM
Clonality	Monoclonal
Purity	Purified by sequential steps of physicochemical fractionation (differential precipitation and solid-phase chromatography methods).
Uniprot ID	075326
Entrez	8482
Dilution Range	Flow cytometry: Recommended dilution: 4 µg/ml. Western blotting: Non-reducing conditions.
Expiration Date	12 months from date of receipt.



Flow
cytometry
analysis
(surface
stainin...