

Product Datasheet

Human MXRA5 ELISA Kit (orb440985)

Description	Human Matrix Remodelling Associated Protein 5 (MXRA5) Sandwich ELISA Kit
Reactivity	Human
Range	0.312-20ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma and other biological fluids.
Sample Types	serum, plasma and other biological fluids.
Uniprot ID	Q9NR99
Sensitivity	0.109ng/mL
Entrez	25878
Expiration Date	Please enquire.