

Product Datasheet

Human MYF5 ELISA Kit (orb440995)

Description	Human Myogenic Factor 5 (MYF5) Sandwich ELISA Kit
Reactivity	Human
Range	0.156-10ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	tissue homogenates and other biological fluids.
Sample Types	tissue homogenates and other biological fluids.
Uniprot ID	P13349
Sensitivity	0.066ng/mL
Entrez	4617
Expiration Date	Please enquire.