

Product Datasheet

Human MYH11 ELISA Kit (orb440997)

Description	Human Myosin Heavy Chain 11, Smooth Muscle (MYH11) Sandwich ELISA Kit
Reactivity	Human
Range	0.781-50ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma and other biological fluids.
Sample Types	serum, plasma and other biological fluids.
Uniprot ID	P35749
Sensitivity	0.27ng/mL
Entrez	4629
Expiration Date	Please enquire.