

Product Datasheet

Mouse MYH1 ELISA Kit (orb440999)

Description	Mouse Myosin Heavy Chain 1, Skeletal Muscle, Adult (MYH1) Sandwich ELISA Kit
Reactivity	Mouse
Range	15.625-1000pg/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma, tissue homogenates and other biological fluids.
Sample Types	serum, plasma, tissue homogenates and other biological fluids.
Uniprot ID	Q5SX40
Sensitivity	7.3pg/mL
Entrez	17879
Expiration Date	Please enquire.