

Product Datasheet

Rat MYO ELISA Kit (orb441029)

Description	Rat Myoglobin (MYO) Sandwich ELISA Kit
Reactivity	Rat
Range	0.625-40ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma, tissue homogenates and other biological fluids.
Sample Types	serum, plasma, tissue homogenates and other biological fluids.
Uniprot ID	Q9QZ76
Sensitivity	0.264ng/mL
Entrez	59108
Expiration Date	Please enquire.