

Product Datasheet

Rabbit MYO ELISA Kit (orb441030)

Description	Rabbit Myoglobin (MYO) Sandwich ELISA Kit
Reactivity	Rabbit
Range	78.125-5000pg/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma and other biological fluids.
Sample Types	serum, plasma and other biological fluids.
Uniprot ID	P02170
Sensitivity	28pg/mL
Entrez	100351487
Expiration Date	Please enquire.