

Product Datasheet

Human NAGLU ELISA Kit (orb441043)

Description	Human N-Acetyl Alpha-D-Glucosaminidase (NAGLU) Sandwich ELISA Kit
Reactivity	Human
Range	0.312-20ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma, tissue homogenates and other biological fluids.
Sample Types	serum, plasma, tissue homogenates and other biological fluids.
Uniprot ID	P54802
Sensitivity	0.125ng/mL
Entrez	4669
Expiration Date	Please enquire.